REF (Catalogue number)	Name of product	Contents
ACR-9902-1	LabStrip U mALB/CREA	25 reagent strips

#### Summary:

25 test strips for the semi-quantitative determination of albumin and creatinine in fresh human urine [1] and the screening of samples for microalbuminuria [2].

Can only be used for evaluation with the help of a DocUReader 2 Pro urine analyser or for visual reading. Under Directive 98/79/EC, LabStrip U mALB/CREA is classified as an in vitro diagnostic medical device intended for professional use.

#### Clinical utility:

The test strips are intended to be used as a screening test for the detection of the pathological changes associated with kidney diseases.

Microalbuminuria is a symptom of the first stage of nephropathy and can also indicate cardiovascular diseases [3], [4]. The concentration of albumin in the urine depends on the amount of urine excreted [5]. The level of creatinine excreted into the urine is generally constant. By using the albumin/creatinine ratio, the volumetric effects can be corrected for; therefore, fresh urine can be used for diagnosing and tracking/monitoring the treatment of microalbuminuria instead of urine collected 24 hours ago [6].

#### Measurement principles:

Each test strip contains 2 test zones, which incorporate sensitive chemicals. If these zones come into contact with urine, a chemical reaction takes place, which results in the discoloration of the test zone.

**Albumin:** The test is based on the so-called 'protein error of indicators' phenomenon, the indicator being a tetrabromophenol-sulfonephthalein derivative in this case [7]. In an acidic environment, the dye binds to the albumin, causing the color of the test strip to change from light to dark turquoise.

**Creatinine:** The test is based on the peroxidase-like activity of a copper-creatinine complex. This complex acts as a catalyst for the color reaction, changing the color of the test zone from light green to dark teal [8].

#### Contents of the package:

The reagent kit suitable for performing 25 tests contains the following items:

- 25 LabStrip U mALB/CREA test strips;
- a label displaying a color scale for visual readings;
- 1 registration card to provide settings for DocUReader 2 Pro device;
- 1 package insert.

#### Other equipment required for the urinalysis process:

- DocUReader 2 Pro urine analyzer device and its User Manual;
- Clean, chemical-free and dry vessel for collecting the urine sample.

#### Reagent composition:

The reagents of the individual test zones contain the following reactive substances:

<b>Albumin:</b>	Tetrabromophenol-sulfonephthalein derivative	1.6 %
<b>Creatinine:</b>	Copper sulphate	1.5 %
	Cumolhydroperoxide	4.0 %
	Tetramethylbenzidine	1.7 %

The reagent composition values are to be interpreted as mass concentration (w/w) specifications used during manufacturing and they may vary within the manufacturing tolerance ranges.

#### Attention!

- All items included in the package can be disposed of as communal waste. As the reactive substances are present in extremely low quantities, they do not fall within the scope of the EU regulations on hazardous materials.
- Ingestion of these substances, as well as their contact with skin or the mucous membranes should be avoided.
- For in vitro diagnostic use only!
- The vial cap containing a non-toxic, molecule-sieve based drying agent which protects the test strips from air-humidity. In the event of accidental ingestion, drink plenty of fluid.
- If you have any questions, please do not hesitate to contact your local distributor.

#### Preparatory steps before urinalysis:

#### Attention!

If you intend to use a medical device for urinalysis, please remember that LabStrip U mALB/CREA strips can only be used with DocUReader 2 Pro urine analyzer devices.

When opening a new package of LabStrip U mALB/CREA strips, you will find 1 registration card, which is required for the settings of the DocUReader 2 Pro device. The information necessary for the strip registration process can be found in the User Manual of the device. The strip registration does not necessary again until the opening of the next strip container.

Please follow the instructions provided in the User Manual of the DocUReader 2 Pro urine analyzer device.

#### Specimen Collection and Preparation:

- Use a clean, chemical-free (free of cleaning agents, disinfectants, or preservatives), and dry vessel for sample collection, which is large enough to allow all reagent zones of the strip to be submerged in the sample.
- Do not add any preservatives to the sample.
- Perform the test on a fresh, thoroughly mixed, uncentrifuged sample. Samples taken during the first morning urine are the most suitable for analysis. It is unnecessary to use urine collected 24 hours earlier [9].
- Protect the sample from direct sunlight or heat.
- Perform the analysis of the sample within 2 hours of collection.
- If immediate analysis is unfeasible, store the sample refrigerated (between +2 and +4°C) but let it warm to room temperature (between +15 and +25 °C) before the start of the analysis process.
- The pH value of the sample may change if stored at room temperature due to the increase in the amount of bacteria in preservative-free urine, which may interfere with the detection of albumin.
- Residues of disinfectants containing chlorhexidine may also interfere with the test results.

#### Steps of the analysis process:

- Remove the required quantity of strips from the container then immediately reseal the container with the lid, which contains a desiccant.
- Do not touch the test zones of the strips.

#### In case of instrumental reading:

- Please read the User Manual of the DocUReader 2 Pro device carefully.
- The DocUReader 2 Pro urine analyzer device detects any changes in color after a 1-minute incubation period and displays the results on the screen.

#### In case of visual reading:

- Submerge the test strip in the thoroughly mixed urine sample for a short period of time (approx. 2 seconds). Ensure that all test zones are submerged in the sample.
- Remove excess urine by touching the side of the strip either to the edge of the vessel (test tube) or to a piece of absorbent paper.
- Hold the test strip horizontally during the incubation period in order to avoid interactions between the test zones located next to each other.
- After 60 seconds from submersion, read the results off the color scale printed on the container. Do not try to read the results if at least 2 minutes have passed since submersion.
- The colors shown on the color scale represent the nominal values of the test zones. The actual values are shown around the nominal values.
- Evaluate the albumin/creatinine ratio based on Table No. 2.
- The compensation test zone is only used in case of instrumental reading.

#### Attention!

- Never remove the strip from the device during the incubation period or the reading process.
- Before the start of the analysis process, please ensure that you have followed the steps included in the User Manual of the DocUReader 2 Pro urine analyzer device precisely.
- DO NOT perform urinalysis if the temperature is below +15°C or above +32°C.
- Never expose the test strips to direct sunlight or heat.
- DO NOT ever reuse a test strip.; Dispose each test strip carefully after use.
- Always store the reagents in their original packaging until they are used; DO NOT mix the strips found in different packages together.
- Never use a test strip for urinalysis if it was left on open air (outside the container) for more than 5 minutes.**
- Never use a test strip past its expiration date.**
- Never use a test strip after 3 months from the first opening of its container.**

#### Biological hazard!

All samples and equipment that may come into contact with urine carry a risk of contamination. Following the completion of the analysis, dispose of them accordingly, in line with the applicable local regulations.

- In certain cases, changes in the test environment or the heterogeneity of urine samples (caused by variations in activators, inhibitors, or ion concentrations) may result in differences in color intensity or contrast.
- The effect of various types of medication, as well as their components and metabolites, on the test is not known in all cases. The color reactions of the test zones may differ. Therefore, in the event of uncertainty, we recommend that you repeat the urinalysis process after the patient has stopped using the relevant medication.
- Always follow local laboratory rules and regulations.
- The test strips DO NOT contain toxic materials.

#### Results and measurement ranges:

When obtaining a visual reading, compare the test strip with the color scale printed on the container label.

**Albumin:** Protein normally is not present in urine.

**Creatinine:** Urine always contains creatinine.

Table No. 1: Measurement ranges of the albumin and creatinine test zones:

Parameters	Categories				
	10	30	80	150	500
Albumin (mg/L)	10	30	80	150	500
Creatinine (mg/dL)	10	50	100	200	300
Compensation pad	Only used in case of instrumental reading.				

Following the determination of the albumin and creatinine concentrations, the albumin/creatinine ratio can be estimated with the help of the interpretation table found on the container label: its value is either 'Normal', 'Abnormal', or 'High Abnormal'.

Table No. 2: Interpreting the albumin/creatinine ratio

		Creatinine				
		A	B	C	D	E
		(10 mg/dL) (0.9 mmol/L)	(50 mg/dL) (4.4 mmol/L)	(100 mg/dL) (8.8 mmol/L)	(200 mg/dL) (17.7 mmol/L)	(300 mg/dL) (26.5 mmol/L)
Albumin	1 (10 mg/L)	X	Normal	Normal	Normal	Normal
	2 (30 mg/L)	High abnormal	Abnormal	Normal	Normal	Normal
	3 (80 mg/L)	High abnormal	Abnormal	Abnormal	Abnormal	Normal
	4 (150 mg/L)	High abnormal	High abnormal	Abnormal	Abnormal	Abnormal
	5 (500 mg/L)	High abnormal	High abnormal	High abnormal	Abnormal	Abnormal

[X]: The sample is too diluted. Please repeat the test using a new sample.

#### Examples for interpreting the results:

- If the urine sample contains 10 mg/L albumin and 100 mg/dL creatinine, the albumin/creatinine ratio is considered normal.
- If the results of the test are as follows: 150 mg/L albumin and 200 mg/dL creatinine, this is considered abnormal (microalbuminuria).

If you are using a DocUReader 2 Pro device, please read its User Manual carefully.

#### Limitations:

Detergents, cleaning agents, disinfectants, and preservatives can cause incorrect albumin or creatinine concentration readings. A high concentration of certain materials typically found in urine, in particular hemoglobin, riboflavin, and bilirubin, may cause atypical discoloration of the test zones.

Never make diagnostic or therapeutic decisions on the basis of results obtained from a single test. When establishing a diagnosis and prescribing the necessary treatment, always evaluate the results obtained from the test strips together with other medical findings related to the patient.

The effect of various types of medication, as well as their components and metabolites on the test is not known in all cases. Therefore, in the event of uncertainty, we recommend that you repeat the urinalysis process after the patient has stopped using the medication in question. Patients may stop using medication only if instructed to do so by their doctor.

Since the composition of urine is not constant (e.g., its activator or inhibitor content, ion concentration, or color may change from sample to sample) and the reaction environment is not always identical, color intensity and shade may differ between samples. If the urine sample is dark in color, this may interfere with the evaluation of the test zones.

#### Interferences:

The following substances do not interfere with the test system up to the concentration specified next to their name: Glucose 1500 mg/dL; Bilirubin 4 mg/dL; Ascorbic acid 80 mg/dL; Urobilinogen 12 mg/dL; Sodium nitrite 1 mg/dL; Leukocytes 500 Leu/μL; Haemoglobin 0.87 mg/dL; Lithium acetoacetate 300 mg/dL; Acetylsalicylic acid 60 mg/dL; Caffeine 60 mg/dL; Riboflavin 10 mg/dL; Creatinine 100 mg/dL; Creatine 10 mg/dL; Formaldehyde 0.1 %; Fructose 100 mg/dL; Galactose 80 mg/dL; Uric acid 50 mg/dL; Urea 400 mg/dL; Lactose 10 mg/dL; Human IgG 5 mg/dL; Myoglobin 0.26 mg/dL. If the substances listed here are found in the urine in a concentration higher than specified, this may interfere with the detection of the albumin/creatinine ratio.

Extremely alkaline samples can cause false positive results to be displayed on the albumin test zones. Other proteins that may be present in the urine may also cause false positive results to be displayed on the albumin test zones.

#### Expected values:

The albumin/creatinine ratio can be specified using mg/g or mg/mmol as the unit of measurement.

Table No. 3: Expected values of the albumin/creatinine ratio:

Classification	Conventional unit (mg/g)	SI unit (mg/mmol)	Interpretation
Normal	≤30	≤3.4	
Abnormal	31-299	3.5-33.8	Microalbuminuria
High Abnormal	≥300	≥33.9	Microalbuminuria, Proteinuria

#### Storage and stability:

The test strips are to be stored in a cool, dry, and dark place in their sealed original packaging at a temperature between +2 and +25°C. After removing a test strip, reseal the container immediately. Do not remove the desiccant from the original cap.

Protect the test strips from direct sunlight, moisture, high temperatures, and chemicals. If stored in accordance with these instructions, the test strips will retain their stability until their expiration date. The test strips retain their stability for 3 months after first opening of the container. Do not touch the reaction surface.

#### System quality control:

Control measurements are recommended to be performed after opening a new test strip container or when beginning to use a new lot of test strips.

During the quality control process, the solutions used, the frequency of the tests, the documentation of the target values, the threshold values, and the results are subject to the internal policies of the testing laboratory as well as the applicable local rules and regulations. If the results fall outside the target range, calibrations must be performed. It is recommended for each laboratory to develop and specify its own quality control policies.

DO NOT use water as a negative control substance.

For regular quality control of the test strips and the devices, we recommend the following urine control solutions: CombiScreen® Dip Check (Analyticon Biotechnologies AG, REF:93010), CombiScreen® Drop Check (Analyticon Biotechnologies AG, REF:93015) or Microbium® Microalbumin Control / Level 1&2 (Quantimetrix Corporation, REF:1250-01). Besides these, the appropriate control solutions of other manufacturers can also be used.

Submerge the strip in the control solution instead of urine. Follow the instructions provided in the Chapter titled 'Steps of the analysis process'.

Please follow the instructions provided in the Chapter of the DocUReader 2 Pro urine analyzer device's User Manual titled 'Quality Control Testing'.

#### Specific performance characteristics:

The performance characteristics of the LabStrip U mALB/CREA urine test strips were established based on clinical and analytical studies. Sensitivity depends on the color detection capability of the device, the presence or absence of any interfering characteristics, and in the case of visual readings, the light conditions. For the purpose of performance evaluation, we used residual urine samples originating from hospital laboratories. We analyzed the samples with LabStrip U mALB/CREA urine test strips, as well as an immunochemical albumin determination technique and the kinetic Jaffe creatinine determination method [10].

#### Accuracy

The test strip typically measured values identical to those obtained with the comparison methods.

Table No. 4 shows the results of the performance evaluation.

Table No. 4: The performance of the LabStrip U mALB/CREA urine test strip in relation to the comparison methods:

Result	a) Accuracy	b) Diagnostic sensitivity (PPA)	c) Diagnostic specificity (NPA)
Albumin/Creatinine ratio (ACR) obtained with visual reading	88 % (n=1261)	92 % (n=1261)	79 % (n=1261)
Albumin/Creatinine ratio (ACR) obtained with instrumenta reading	85 % (n=1264)	93 % (n=1264)	70 % (n=1624)

#### Precision

Intra assay (within-run) precision (repeatability) was established by performing 20 consecutive measurements on each of two control urine samples (one normal and one abnormal), using 3 different lots of LabStrip U mALB/CREA urine test strips.

Inter assay (day-to-day) precision (reproducibility) was established by performing 20 measurements on each of two control urine samples (one normal and one abnormal) on 20 different days. Over the course of the 20-day long experiment, 3 different lots of reagents were used.

Table No. 5: The precision of LabStrip U mALB/CREA urine test strips:

Visual readings:					
Albumin	Intra assay	Inter assay	Creatinine	Intra assay	Inter assay
10 mg/L	100 %	100 %	50 mg/dL	95 %	90 %
80 mg/L	95 %	90 %	200 mg/dL	95 %	90 %

#### Using the device:

Albumin	Intra assay	Inter assay	Creatinine	Intra assay	Inter assay
10 mg/L	100 %	100 %	50 mg/dL	95 %	90 %
80 mg/L	95 %	90 %	200 mg/dL	95 %	90 %

#### Analytical sensitivity

**Visual readings:** The test indicates the presence of albumin or creatinine if the sample contains at least 20 mg/L of albumin or 10 mg/dL of creatinine.

**Using the device:** The test indicates the presence of albumin or creatinine if the sample contains at least 20 mg/L of albumin or 10 mg/dL of creatinine.

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#### Bibliography / Bestellnummer / Bibliografía:

- [1] Parsons M., Newman D.J., Pugia M., Newall R.G., and Price C.P.: *Performance of a reagent strip device for quantitation of the urine albumin : creatinine ratio in a point of care setting.* Clin Nephrol. 1999 Apr;51(4):220-7.
- [2] Sasaki M., Pugia M.J., Parker D.R., Kuromoto K., Furukawa I., and Konishi I.: *Measurement of the albumin content of urinary protein using dipsticks.* J Clin Lab Anal. 1999;13(5):246-50.
- [3] Kaplan N.M.: *Microalbuminuria: a risk factor for vascular and renal complications of hypertension.* Am J Med. 1992 Apr 27;92(4B):85-125.
- [4] Apostolovic S., Stanojevic D., Djordjevic V., Jankovic Tomasevic, R., Salinger Martinovic S., Djordjevic Radojkovic D., Koracevic G., Stojanovic I., and Pavlovic M.: *Prognostic significance of microalbuminuria in patients with acute myocardial infarction.* Clin. Lab. 57; 229-235; 2011
- [5] Pugia M.J., Lott J.A., Clark L.W., Parker D.R., Wallace J.F., Willis T.W.: *Comparison of urine dipsticks with quantitative methods for microalbuminuria.* Eur J Clin Chem Clin Biochem. 1997 Sep;35(9):693-700.
- [6] **Position Statement: Nephropathy in Diabetes.** Diabetes Care, 27, S79-S83, 2004
- [7] Pugia M.J., Lott J.A., Profitt J.A., and Cast T.K.: *High-sensitivity dye binding assay for albumin in urine.* J Clin Lab Anal. 1999;13(4):180-7.
- [8] Pugia M.J., Lott J.A., Wallace J.F., Cast T.K., Bierbaum L.D.: *Assay of creatinine using the peroxidase activity of copper-creatinine complexes.* Clin Biochem. 2000 Feb;33(1):63-70.
- [9] Guy M., Borzomato J.K., Newall R.G., Kalra P.A., and Price C.P.: *Protein and albumin-to-creatinine ratios in random urines accurately predict 24 h protein and albumin loss in patients with kidney disease.* Ann Clin Biochem. 2009;46: 468-76.
- [10] Burtis C.A., Ashwood E.R. and Bruns D.E. (eds): *Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (5th edition).* Elsevier, St. Louis, USA, 2012, 2238 pp, 909 illustrations. ISBN: 978-1-4160-6164-9.

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